



SpinStar™ Total DNA Kit

Extraction of Total DNA from a
Broad Range of Specimens

Extraction from a Wide Range of Samples

The SpinStar™ Total DNA Kit is intended for purification of total DNA (e.g. genomic, mitochondrial and of bacterial, viral and parasite origin) from up to 400 µl human whole blood, serum, plasma, respiratory specimen in VTM and up to 200 µl sputum (pretreated) or stool suspension.

The kit offers efficient, fast and convenient workflows for the **broad range of validated sample materials** and consistent high-quality DNA through proven silica membrane technology in spin column format.

The purified DNA is ready to use for downstream analysis like PCR, real-time PCR or any other molecular application that requires pure DNA of highest quality and integrity.



Order Information

Product	Content	Σ	Cat. No.
SpinStar™ Total DNA Kit 2.0	SpinStar™ columns, collection tubes, lysis buffer, genomic DNA binding buffer, wash buffer 1, wash buffer 2, elution buffer, lysis enhancer, Proteinase K, user manual	100	821803

AstronDX Technologies

Unit 307 Block B, Phileo Damansara 1, 9 Jalan 16/11,
46350 Petaling Jaya, Selangor, Malaysia

phone +603 7931 6760
fax +603 7931 5352
email info@astrondx.com
website www.astrondx.com

