




# SpinStar™ Viral Nucleic Acid Kit

Simultaneous Extraction of Viral  
Nucleic Acids (DNA and RNA)

## Efficient Extraction of Viral Nucleic Acids

The SpinStar™ Viral Nucleic Acid Kit is intended for simultaneous purification of viral nucleic acids from fresh or frozen human serum and plasma, respiratory specimens in VTM and stool.

The kit is a validated one stop solution for extraction of **both viral DNA and viral RNA in a single workflow** and offers fast, convenient and consistent high-quality results through proven silica membrane technology in spin column format. The purified DNA and/or RNA is ready to use for downstream analysis like PCR, real-time PCR or any other molecular application that requires pure nucleic acids of high quality and integrity.



## Order Information

Product	Content	Σ	Cat. No.
SpinStar™ Viral Nucleic Acid Kit 1.0	SpinStar™ columns, collection tubes, lysis buffer, wash buffer 1, wash buffer 2, elution buffer, carrier RNA, Proteinase K, user manual	100	811803

## AstronDX Technologies

Unit 307 Block B, Phileo Damansara 1, 9 Jalan 16/11,  
46350 Petaling Jaya, Selangor, Malaysia

phone +603 7931 6760  
fax +603 7931 5352  
email [info@astrondx.com](mailto:info@astrondx.com)  
website [www.astrondx.com](http://www.astrondx.com)

