



# SpinStar™ VNAplus Mix

For Use in Combination with the  
SpinStar™ Viral Nucleic Acid Kit

## Ultra-Sensitive Extraction of Viral Targets

The SpinStar™ VNAplus Mix allows concentration of viral DNA and RNA from 1 ml plasma or serum, prior to purification with the SpinStar™ Viral Nucleic Acid Kit. Ultracentrifugation and specialized laboratory equipment are not required.

The SpinStar™ VNAplus Mix forms complexes with nucleic acids, which can be sedimented by centrifugation to form a pellet. This pellet is then resuspended in standard volume of lysis buffer and nucleic acids are subsequently purified using the SpinStar™ Viral Nucleic Acid procedure, enabling efficient and ultra-sensitive viral nucleic acid extraction.



## Order Information

Product	Content	Σ	Cat. No.
SpinStar™ VNAplus Mix	90 ml ready-to-use VNAplus solution	100	812803

## AstronDX Technologies

Unit 307 Block B, Phileo Damansara 1, 9 Jalan 16/11,  
46350 Petaling Jaya, Selangor, Malaysia

phone +603 7931 6760  
fax +603 7931 5352  
email [info@astrondx.com](mailto:info@astrondx.com)  
website [www.astrondx.com](http://www.astrondx.com)

